



Product Code . JLI-CLE-10553

Forensics Using Simulated Blood

Description

Though the use of blood type in a forensic investigation is not enough to prove guilt, JLab India it may aid in exonerating a potential suspect. In this activity, students act as lab technicians and assist investigators by examining evidence collected at a crime scene. Students first use a

presumptive blood test to determine if a substance on a stained piece of cloth may be blood, and then determine the blood type of a sample collected at the scene. Students then compare their results to those of samples provided from two suspects. Kit contains enough materials for 10 groups. Teacher's Manual and Student Study Guide copymasters are included.

Kit Includes:

1 btl Blood detection reagent #1

1 btl Blood detection reagent #2

4 Simulated blood samples

40 Blood typing trays

10 Cotton swabs

1 pc Cloth

Suspect #1

Suspect #2

Crime Scene Evidence




1 btl Simulated anti-A serum

1 btl Simulated anti-B serum

1 btl 70% ethanol

Victim

1 btl Simulated anti-Rh serum

Educational Lab Equipments,
#449, HSIIDC, Industrial Area, Saha, Haryana
Direct Contact Details  +91-98173-19615  sales@educational-equipments.com
 www.educational-equipments.com