



Product Code . JLI-BLE-10575

## Tulgren Funnel

### Description

---

Description :-

A soil sample is placed on the perforated plate and a 220V heating lamp is held above the soil, driving the organisms down through the funnel into a collecting jar. Supplied with coarse and fine mesh. Horizontal support rod and bosshead included.

Technical Specifications :-

Sieve diameter 170mm




Sieve depth 65mm

Funnel diameter 175mm

---

Requires a low wattage bulb (B22D, Bayonet fitting 22mm diameter, 40 watt or less).

---

Educational Lab Equipments,  
#449, HSIIDC, Industrial Area, Saha, Haryana  
Direct Contact Details  +91-98173-19615  [sales@educational-equipments.com](mailto:sales@educational-equipments.com)  
 [www.educational-equipments.com](http://www.educational-equipments.com)