



Educational Lab Equipments



Description


Examine some possible methods of gathering evidence when dusting for fingerprints is not effective. Learn to identify fingerprint types, a method of fuming for fingerprints, and a technique of chemically-developing fingerprints.


Utilise alternative methods for detecting fingerprints.


The three activities include fingerprint analysis, ninhydrin development, and cyanoacrylate fuming.

This activity includes a Teachers Manual and Student Guide and Analysis copy-masters. There is enough material for 15 groups.

Educational Lab Equipments,
#449, HSIIDC, Industrial Area, Saha,
Haryana

Direct Contact Details 

+91-98173-19615  sales@educational-equipments.com

 www.educational-equipments.com